

## How to Use the EpiPen Auto-Injector

# Blue to the sky. Orange to the thigh.

---

### How to use EpiPen® and EpiPen® Jr (epinephrine) Auto-injectors.

Remove the EpiPen® Auto-Injector from the carrier tube and follow these 2 simple steps:



- Hold firmly with orange tip pointing downward.
- Remove blue safety cap by pulling straight up. Do not bend or twist.



- Swing and push orange tip firmly into mid-outer thigh until you hear a 'click'.
- Hold on thigh for several seconds.



- Built-in needle protection**
- After injection, the orange cover automatically extends to ensure the needle is never exposed.